



**Sequence polymer flanking the insert - Multiple cloning site removed**

AATTGGAGGCTACAGTCAGTGGAGAGGACTTTC \*\*\*gene\*\*\* GAAAGAGGAAGGGCTGGAAGAGGAAGGAGCT

\* \* \* \* \*

3141 3171 1 33

\*\*\* Your gene may be flanked by a unique sequence of up to 18 bp on each side. If your gene is cloned with a flanking sequence an information sheet with the sequence will be included with your gene shipment and electronically posted on your GeneMaker Secure Website File Area in the "Messages from Blue Heron" folder. We recommend ordering your sequence with flanking restriction sites to facilitate future subcloning.